
Vitrification Cryotop Method for Oocyte and Embryo

CE 2797

Code: /T801

Symbols on the Label

Ver.13

	Sterilized using aseptic processing techniques
	Do not reuse
	Catalog number
	Batch code
	Use by
	Storage temperature
	CE mark
	Manufacturer
	Do not re-sterilize
	Do not use if package is damaged
	Consult instructions for use
	Authorized representative in the European Community

EN: Intended use
This product is to be used for vitrification of oocytes (MI) and embryos.

Vitrification Media

- No. 0 (Color code: white) Basic Solution (BS): 1x1.5 ml vial (only for oocyte (MI) vitrification)
- No. 1 (Color code: blue) Equilibration Solution (ES): 1x1.5 ml vial
- No. 2 (Color code: green) Vitrification Solution (VS): 2x1.5 ml vials

⚠ **Caution:** Before use, check the specifications of the container and labeling (number marked on top of the cap, cap color, vial labeling color, solution name, and volume). If you notice anything unusual regarding the items mentioned above, do not use the product and please contact the distributor.

Recommended equipment

- Cryotop: 1 Cryotop stores up to 4 oocytes (MI) or 4 Embryos as a recommendation.
- ReproPlate with 6 wells

⚠ **Caution:** This product is sterilized. Please handle the product in a sterile field with a clean bench.

Instructions for use (IFU)

Preparation

- Fill 90% of the Cooling Rack with fresh liquid nitrogen. Compare the width of petriwell space with the thickness of zona pelliculosa and record it.
- Use a sterilized micropipette as a handling tool, with a suitable internal diameter for oocytes or embryos. The recommended internal diameters are as follows: 120 µm for oocytes (MI), 120 µm for pronuclear stage embryos, 150 to 180 µm for cleavage stage embryos, and 180 to 250 µm for blastocyst stage embryos.

Equilibration

⚠ **Caution:** Equilibration procedures for Oocyte (MI) and Embryo are different.

Equilibration of oocytes (MI)

- Drop 20 µl of BS and 300 µl each for VS1 and VS2 on the Repro Plate with micropipette.
- Transfer the Oocyte (MI) from the culture dish to the BOTTOM of BS.
- Step 1. Add ES 20 µl gently on the TOP of BS and leave it for 3 minutes.
- Step 2. Add another ES 20 µl gently on the TOP of BS and leave it for 3 minutes.
- Step 3. Add another ES 240 µl gently on the TOP of BS and leave it for 9 minutes.

Equilibration of embryos

- Drop 300 µl for ES, VS1 and VS2 on the Repro plate using micropipette.
- Transfer the Embryo to the TOP center of ES.
- It will spontaneously begin to swim and then gradually return to its original size by absorbing the ES solution (within 15 minutes).

Vitrification

⚠ **Caution:** The following steps from 1 to 9 should be completed between 60 and 90 seconds.

- Aspirate the oocyte (MI) embryo from ES with the tip of a micropipette.
- Transfer the oocyte (MI) embryo to the TOP center of VS1.
- Aspirate the oocyte (MI) embryo with micropipette and blow it out. Repeat this process 3 times, changing the position in VS1.
- Transfer the oocyte (MI) embryo to VS2.
- Change the position of the oocyte (MI) embryo in VS2 with a micropipette.
- Place the oocyte (MI) embryo by the black line on the Cryotop.
- Make a planar droplet.
- Apply the oocyte (MI) embryo on the Cryotop with a minimal volume of VS2 (less than 0.1 µl) under a microscope.
- Plunge the Cryotop immediately in liquid nitrogen.
- Put the Cryotop in a can and store it in a storage tank.

Quality control specification

The following tests were performed for each of this product:

- Sterility by the current USP Sterility Test <71> (Solutions)
- Endotoxin by LAL methodology
- Mouse Embryo Assay (MEC)
- pH (Ph. Eur., USP)
- Osmolality (Ph. Eur., USP)

Storage instructions and stability

Solutions: Store the vials at 2 to 8 °C.
This product is stable until the expiry date labeled on the vial.

Composition

- HEPES within Basic Culture Medium
- Ethylene Glycol
- Dimethyl Sulfoxide
- Trehalose
- Hydroxypropyl Cellulose

Contraindications

- Do not re-sterilize.
- Do not re-use.
- Do not use solution that shows cloudiness or has turned yellow.
- Wear sterile gloves, hat, white coat to prevent contamination
- Use the product if you notice anything unusual regarding the specifications on the label (number, color, name, volume).
- Sterile if the package or container is unopened or undamaged. Do not use if the package or container are opened or damaged.
- Upon delivery media must be stored in original unopened container and refrigerated at 2 to 8 °C.
- Do not use and please discard if the media is not stored under refrigeration (2 to 8 °C).
- Do not use the product if past the expiration date.

Warning

Read the instructions for use prior to use.

Use the thawing media T802 to thaw the oocytes (MI) or embryos vitrified with KITAZATO vitrification media

This product is intended to be used by medical specialists trained in fertility treatment.

Aseptic technique should be used.

Use sterilized equipment and materials only.

Decontaminate the workroom for aseptic procedures.

Follow procedures in an environmentally controlled room.

Confirm the level of liquid nitrogen once every two weeks.

Wear sterilized gloves, hat, white coat to prevent contamination

Operate must be conducted by two people and perform a double check

In case of eye or skin contact with Vitrification/Thawing media, immediately flush eye/skin with water.

Morphologically abnormal oocytes, embryos, or significantly poor grade oocytes or embryos are unsuitable for cryopreservation

Operation must be conducted by two people and perform a double check

Observe all federal, state and local environmental regulations when handling the product.

Request a qualified waste disposal contractor to dispose the products as residual and non-reusable solutions.

Dispose the contaminated containers and packages the same way as instructed for containers containing the product inside.

In case of infection, dispose the product appropriately in a prescribed manner.

The user shall be responsible for any problems caused by non-conformity to the present IFU.

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Copyright © Kitazato Corporation All Rights Reserved.