CRYOTOP KITAZATO

Copyright © KITAZATO CORPORATION. All Rights Reserved.

Kitazato Corporation 100-10 Yanagishima, Fuiji, Shizuoka 416-0932, Japan

TEL: +81-545-65-7122 FAX: +81-545-65-7128

Dibimed - Biomedical Supply S.L. C/Jorge Corin, 3 Bajo 1-2, 46015 Valencia, Spain

TEL: +34-96-305-63-95 FAX: +34-96-305-63-96

CE 2797

ENI CRYOTOP

INTENDED USE: Cryotopis used for storage of vitrified human oocytes and embryos.

QUALITY CONTROL TESTING: Each lot of Cryotop receives the following tests:
• Sterility by the current USP Sterility Test
• Endotoxin by LAL methodology
• Mouse Embryo Assay (one cell).
All results of each lot are reported on a Certificate of Analysis which is available upon request.

STORAGE INSTRUCTIONS AND STABILITY: Store in original sterile pack at 15-30°C. Cryotop is stable until the expiration date shown on package label when stored as directed.

PRECAUTIONS AND WARNINGS: Sterile if the package is unopened or undamaged. Do not use if package is broken. Each Cryotop is intended for single use. Perform all vitrification and thawing operations for oocytes or embryos at room temperature (23-27 °C). Liquid nitrogen may cause freeze burns if in contact with skin or eyes. Use appropriate precautions when work with liquid nitrogen. Do not re-use. Re-use may cause change in product quality and raise the risk of poor embryo survival rate. This product is intended to be used by medical specialists trained in fertility treatment. Aseptic technique should be used.

INSTRUCTIONS FOR USE:

VITRIFICATION 1. Open the sterile pack of Cryotop and write necessary information about a patient on the handle of Cryotop. 2. Fill 90% of a container with fresh liquid nitrogen. 3. Prepare embryos or oocytes for vitrification according to laboratory protocol. 4. Take the straw cap of Cryotop and plunge it into the fresh liquid nitrogen. 5. Place the Cryotop under the microscope and adjust the focus on the black mark of the Cryotop. Logo should be up. 6. Gently place the oocytes or embryos by the black part of Cryotop sheet with minimal volume of vitrification solution using a suitable pipette. Notes: The recommended load of the Cryotop is up to 4 oocytes or embryos. 7. Quickly plunge the Cryotop into fresh liquid nitrogen. 8. Hold the straw cap with tweezers and inset the Cryotop. Fit the Cryotop with the cap by hands with screwing tightly. 9. Put the Cryotop in a cane and store it in a tank for long term storage. Caution: Take care that the Cryotop remains immersed in liquid nitrogen at all times until thawing.

THAWING: 1. Prepare the thawing media according to laboratory protocol. 2. Retrieve the cane which has the specific Cryotop and quickly immerse the cane in a container filled with fresh liquid nitrogen. Collect the specific Cryotop from the cane in the liquid nitrogen. Check the information about the patient on the handle of the Cryotop. Caution: The Cryotop except the handle remain immersed in liquid nitrogen at all time. Place a container for liquid nitrogen by the stereo microscope. 3. Carefully twist and remove the straw cap from the Cryotop in liquid nitrogen. 4. Quickly immerse the Cryotop sheet into the thawing solution. It should be one second. 5. One minute after immersing into the thawing solution, gently aspirate the oocyte/embryo using a suitable pipette. 6. Perform the thawing procedure according to laboratory protocol.

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PRECAUTIONS AND WARNINGS: Sterile if the package is unopened or undamaged. Do not use if package is broken. Each Cryotop is intended for single use. Perform all vitrification and thawing operations for oocytes or embryos at room temperature (23-27 °C). Liquid nitrogen may cause freeze burns if in contact with skin or eyes. Use appropriate precautions when work with liquid nitrogen. Do not re-use. Re-use may cause change in product quality and raise the risk of poor embryo survival rate. This product is intended to be used by medical specialists trained in fertility treatment. Aseptic technique should be used.
INSTRUCTIONS FOR USE:
VITRIFICATION 1. Open the sterile pack of Cryotop and write necessary information about a patient on the handle of Cryotop. 2. Fill 90% of a container with fresh liquid nitrogen. 3. Prepare embryos or oocytes for vitrification according to laboratory protocol. 4. Take the straw cap of Cryotop and plunge it into the fresh liquid nitrogen. 5. Place the Cryotop under the microscope and adjust the focus on the black mark of the Cryotop. Logo should be up. 6. Gently place the oocytes or embryos by the black part of Cryotop sheet with minimal volume of vitrification solution using a suitable pipette. Notes: The recommended load of the Cryotop is up to 4 oocytes or embryos. 7. Quickly plunge the Cryotop into fresh liquid nitrogen. 8. Hold the straw cap with tweezers and inset the Cryotop. Fit the Cryotop with the cap by hands with screwing tightly. 9. Put the Cryotop in a cane and store it in a tank for long term storage. Caution: Take care that the Cryotop remains immersed in liquid nitrogen at all times until thawing.
THAWING: 1. Prepare the thawing media according to laboratory protocol. 2. Retrieve the cane which has the specific Cryotop and quickly immerse the cane in a container filled with fresh liquid nitrogen. Collect the specific Cryotop from the cane in the liquid nitrogen. Check the information about the patient on the handle of the Cryotop. Caution: The Cryotop except the handle remain immersed in liquid nitrogen at all time. Place a container for liquid nitrogen by the stereo microscope. 3. Carefully twist and remove the straw cap from the Cryotop in liquid nitrogen. 4. Quickly immerse the Cryotop sheet into the thawing solution. It should be one second. 5. One minute after immersing into the thawing solution, gently aspirate the oocyte/embryo using a suitable pipette. 6. Perform the thawing procedure according to laboratory protocol.

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PRECAUTIONS AND WARNINGS: Sterile if the package is unopened or undamaged. Do not use if package is broken. Each Cryotop is intended for single use. Perform all vitrification and thawing operations for oocytes or embryos at room temperature (23-27 °C). Liquid nitrogen may cause freeze burns if in contact with skin or les yeux. Use appropriate precautions when working with liquid nitrogen. Do not re-use. Re-use may cause change in product quality and raise the risk of poor embryo survival rate. This product is intended to be used by medical specialists trained in fertility treatment. Aseptic technique should be used.
INSTRUCTIONS FOR USE:
VITRIFICATION 1. Open the sterile pack of Cryotop and write necessary information about a patient on the handle of Cryotop. 2. Fill 90% of a container with fresh liquid nitrogen. 3. Prepare embryos or oocytes for vitrification according to laboratory protocol. 4. Take the straw cap of Cryotop and plunge it into the fresh liquid nitrogen. 5. Place the Cryotop under the microscope and adjust the focus on the black mark of the Cryotop. Logo should be up. 6. Gently place the oocytes or embryos by the black part of Cryotop sheet with minimal volume of vitrification solution using a suitable pipette. Notes: The recommended load of the Cryotop is up to 4 oocytes or embryos. 7. Quickly plunge the Cryotop into fresh liquid nitrogen. 8. Hold the straw cap with tweezers and inset the Cryotop. Fit the Cryotop with the cap by hands with screw