

CRYOTOP® METHOD

VITRIFICATION PROTOCOL

Vitrification Media VT601 / VT602

Cryotop® **Open System**

Cryotop® **SC Closed System**



KITAZATO®

VITRIFICATION

—————

VITRIFICATION

Vitrification procedure



It is different procedures between Oocyte and Embryo for Equilibration.

PART 1

MATERIALS REQUIRED

Vitrification VT601

No.0 Basic Solution (BS): 1 X 1.5mL vial (Only for Oocyte Vitrification)

No.1 Equilibration Solution (ES): 1 X 1.5mL vial

No.2 Vitrification Solution (VS): 2 X 1.5mL vials

Cryotop®

Cryotop®

Cryotop®SC

- Repro Plate - K6
- Cooling Rack, Steel, Lid
- Pasteur Pipette ****refer to CAUTION**
- Stereomicroscope (Turn off the heating plate)
- Stopwatch or Timer (with count up function)
- Liquid Nitrogen
- Tweezers
- 2 Micro pipettes: 2-20uL / 100-1000uL
- Cane
- Storage tank

Additional Materials for Cryotop®SC

- Cooling Rack SC
- Straw Cutter
- Aluminum Block
- Sealer

CAUTION

Use a pasteur pipette that has a suitable internal diameter for Oocyte or Embryo. The external diameter of Oocyte is about 120um and for Embryo, about 120-250um. This is to optimize the volume of the solutions for the best dilution condition to get the highest survival rate.



PART 2

PREPARATION FOR VITRIFICATION

1

Bring BS, ES and VS to room temperature (25-27°C).



Figure 2-1

2

Write necessary information about a patient on the handle/straw cap of Cryotop (See Figure 2-1). You can also label them.




Figure 2-2

3

[Cryotop]

Fill 90% of Cooling Rack with fresh liquid nitrogen.

[Cryotop SC for closed system]

Place Aluminum Block in Cooling Rack SC from the beginning. Then fill with fresh liquid nitrogen until it covers the top of the Aluminum Block (See Figure 2-3).





Figure 2-3

4

Remove the culture dish containing Oocyte or Embryo from the incubator. Check the quality of the Oocyte or the Embryo well with pasteur pipette under the microscope (See Figure 2-4).

For Oocyte Vitrification, take the cumulus cells off.

Oocyte



※ Compare the width of perivitelline space with thickness of zona pellucida and record it (Ex. 1:1). It makes easy to know the completing of the equilibration after immersing in ES.


Figure 2-4

VITRIFICATION

PART 3

EQUILIBRATION

Oocyte equilibration



Oocyte Equilibration 1

Write BS, VS1 and VS2 on the lid of Repro Plate. Drop 20uL for BS and 300uL each for VS1 and VS2 on the plate with micro pipette (See Figure 3-1). Immediately put the lid on the Repro Plate.

Oocyte Equilibration 2

Aspirate the Oocyte at the tip of the pasteur pipette. Transfer the Oocyte with minimal volume of medium from the culture dish to the **BOTTOM** of BS (20uL).

Oocyte Equilibration 3 - For 3 minutes

Set up the stop watch (with count up function). Check the time with the stop watch for the following steps. Add ES 20 L gently to the **TOP** of BS with the Oocyte moving micro pipette along the well and leave it for 3 minutes (See Figure 3-2).

Oocyte Equilibration 4 - For 3 minutes

Add another ES 20 L gently to the **TOP** of BS as well and leave it for 3 minutes (See Figure 3-2).

Oocyte Equilibration 5 - For 6 - 9 minutes

Add another ES 240 L gently to the **TOP** of BS and leave it for 6 - 9 minutes (See Figure 3-2).

For Equilibration, the volume of Oocyte is required to be recovered completely. Oocyte Equilibration is complete when the width of perivitelline space becomes equal to the width before immersing to ES.



Figure 3-1





Figure 3-2

Embryo equilibration



Recommendation

2PN , 4-cell or 8-cell: 10 – 12min
 Morula or Blastocyst: 12 – 15min

Embryo Equilibration 1

Write **ES**, **VS1** and **VS2** on the lid of Repro Plate. Gently invert each vial of **ES** and **VS** twice to mix contents. Drop each solution 300uL on the plate using micro pipette (See Figure 3-3). Immediately put the lid on the Repro Plate.

Embryo Equilibration 2

Aspirate the Embryo at the tip of the pasteur pipette (See Figure 3-4). Put the Embryo with minimal volume of medium to the **TOP** center of **ES**.

Embryo Equilibration 3 - For 10 - 15 minutes

Set up the stop watch (with count up function). Check the time with the stop watch for the following steps. The Embryo free-falls within 30 seconds.

It spontaneously begins to shrink and then gradually returns to its original size with infiltrating **ES**, which indicates that the Equilibration is complete.




Figure 3-3




Figure 3-4

CAUTION

For Blastocyst Equilibration, wait for disappearing of the perivitelline space. Especially, for vitrification of Blastocyst, Day 5 is recommended.

Equilibration time is as follows:

Oocyte: 12-15min
 2PN, 4-cell or 8-cell: 10-15min
 Morula or Blastocyst: 12-15min



VITRIFICATION

PART 4

VITRIFICATION

It is the same procedure for Oocyte and Embryo.

Vitrification 1

After the completion of Equilibration, aspirate the Oocyte (Embryo) in ES at the tip of pasteur pipette (See Figure 4-1). Transfer the Oocyte (Embryo) to the surface center of VS1 with minimal volume of ES. Blow only the Oocyte (Embryo) out to VS1. To avoid getting the remaining ES in the pasteur pipette into the VS1, blow out the ES to the outside of the well. Aspirate fresh VS1 and blow it out again to the outside of the well. Aspirate fresh VS1 into the pasteur pipette. (See Figure 4-3).




Figure 4-1

Vitrification 2 - Within 0.5 minute

Aspirate the Oocyte (Embryo) in VS1 with the pasteur pipette and blow it out to VS1. Quickly stir five times around the Oocyte (Embryo). Repeat the aspirating, blowing out and stirring three times around changing the positions in VS1 (See Figure 4-2). Displace the outer solution of the Oocyte (Embryo) to VS1 completely until the remaining ES visually disappears.

Vitrification 3 - Within 0.5 minute

Blow out the remaining VS1 in the pasteur pipette to the outside of the well. Aspirate fresh VS2 into the pasteur pipette, and then aspirate the Oocyte (Embryo) in VS1 at the tip of the pipette. Transfer the Oocyte (Embryo) to VS2 with minimal volume of VS1. Stir around the Oocyte (Embryo) changing positions twice with the pasteur pipette in VS2 (See Figure 4-2). This step is completed when the outer Oocyte (Embryo) is displaced to VS perfectly and the flat shrinking in cause of dehydration is observed.




Figure 4-2

Vitrification 4

Place the Cryotop under a microscope (Logo should be up) and adjust the focus on the black mark of the Cryotop sheet (See Figure 4-3).



Figure 4-3

Vitrification 5

Aspirate the shrunk Oocyte (Embryo) in VS2 at the tip of the pasteur pipette (See Figure 4-4). Place the Oocyte (Embryo) by the black mark of Cryotop sheet with minimal volume (less than 0.1/L) of VS2 (See Figure 4-5a and 4-5b). For more than 2 Oocytes (Embryos), make 1 droplet for each (See Figure 4-6a and 4-6b).





Figure 4-4

Removal of the excess VS on the sheet


After putting Oocytes (Embryos) on the Cryotop sheet, the excess VS should be removed by aspirating using pipette.

Step 1




Put the top of the pipette on the bottom end of the big VS drop.

Step 2




Slide the pipette horizontally to outside, and make the VS drop lower.

Step 3



Aspirate the excess VS, and minimize the VS drop (not aspirating oocyte).


Good example



Make a planar droplet by the black mark of Cryotop sheet.

Figure 4-5a

Bad example



Make a steric droplet by the black mark of Cryotop sheet. The volume of VS2 is too much.

Figure 4-5b

Good example




Figure 4-6a

Bad example





Figure 4-6b

VITRIFICATION


It is the same procedure for Oocyte and Embryo.

Cryotop® - Open System Vitrification 6-A


Punge the Cryotop directly into liquid nitrogen. Hold the straw gap wot wweezers and insert the Cryotop from sheet end Pursid nitrogen. Then f the Cryotop with the straw cap by hands screwing tightly in the air (See Figure 4-7).



Hold the straw cap with tweezers and insert the Cryotop into it.



Hold the straw cap with fingers and fit it.



Twist it and make sure if the straw cap fits tightly to the Cryotop.

Figure 4-7



CAUTION

Keep the Cryotop sheet in the liquid nitrogen until transferring to a storage tank. In transferring the Cryotop to other storage tank, keep it in liquid nitrogen. Do not expose of the Cryotop in air until Thawing.

Cryotop®SC - Complete Glosed System Vitrification 6-B

Without direct contact with liquid nitrogen, insert the Cryotop®SC into the straw cap pre-set at the Aluminum Block. Then seal the straw cap (See Figure 4-8).





Figure 4-8

POINT

Lean the upper part of the straw cap against the Cooling Rack. This positioning avoids influence of cool air from liquid nitrogen.



CAUTION

Keep the CryotopSC in the liquid nitrogen until transferring to a storage tank. In transferring the CryotopSC to other storage tank, keep it in liquid nitrogen. Do not expose of the CryotopSC in air until Thawing.



A series of horizontal blue lines for writing notes, consisting of 20 lines spaced evenly down the page.

THAWING

THAWING

PART 1

MATERIALS REQUIRED

Thawing Media VT602

No. 1 Thawing Solution (TS)

No.2 Diluent Solution (DS)

No.3 Washing Solution (WS)

- Repro Plate - K6
- 2 Petri Dish
- Cooling Rack, Lid, Steelbox
- Pasteur Pipette **refer to CAUTION
- Stereomicroscope (Turn off the heating plate)
- Stopwatch or Timer (with count up function)
- Liquid Nitrogen
- Tweezers
- 1 Micro pipette: 100-1000uL

Additional Materials for Cryotop@SC

- Straw Cutter

CAUTION

Use a pasteur pipette that has a suitable internal diameter for Oocyte or Embryo. The external diameter of Oocyte is about 120um and for Embryo, about 120-250um. This is to optimize the volume of the solutions for the best dilution condition to get the highest survival rate.



PART 2

PREPARATION FOR THAWING

1

Warm **TS** vial (sealed) with a Petri Dish in an incubator to 37°C(>1.5hours).

2

Bring **DS** and **WS** to room temperature (25~27°C).

3

Retrieve the cane which has the specific Cryotop, quickly immerse the cane in a Cooling Rack filled with fresh liquid nitrogen. Retrieve the specific Cryotop from the cane in the liquid nitrogen. Check the information about the patient on the label of Cryotop.



CAUTION

Place the Cooling Rack by the stereo microscope

4

Write **DS**, **WS1** and **WS2** on the lid of a Repro Plate. Gently invert each vial of **DS** and **WS** twice to mix contents. Drop 300µL each for **DS**, **WS1** and **WS2** on the Repro Plate with micro pipette. Place it on the microscope stage and lid it. Remove **TS** vial and the Petri Dish from the incubator and place the Petri Dish on the microscope stage. Gently invert the vial of **TS** twice to mix contents and pour the full contents into the Petri Dish (See Figure 2-1).



Figure 2-1

5

Adjust the focus of the microscope to the Petri Dish with **TS**.
Use pasteur pipette in order to focus easily on the center of the Petri Dish (See Figure 2-2).




Figure 2-2


THAWING

PART 3

THAWING

Cryotop® and Cryotop®SC have different procedures.

Cryotop® - Open System



Thawing 1

Carefully twist and remove the straw cap from the Cryotop in liquid nitrogen (See Figure 3-1). Prop it against the corner of the Cooling Rack.




Figure 3-1

Thawing 2

Be ready to use pasteur pipette keeping the Cryotop in liquid nitrogen. Set up the stop watch (with count up function). Check the time with the stop watch for the following steps.

Set up the stop watch (with count up function). Check the time with the stop watch for the following steps.

Thawing 3 - For 1 minute

Quickly immerse Cryotop sheet into TS on the microscope stage. It should be within 1 second (See figure 3-2). Find the Oocyte (Embryo) adjusting the focus on the black mark of the Cryotop sheet. 1 minute after immersing into TS, gently aspirate the Oocyte (Embryo) with the pasteur pipette after dispensing it from the sheet. Aspirate the Oocyte (Embryo) even if it does not dispense from the sheet. Also, aspirate TS until the Oocyte (Embryo) reaches 2mm from the tip of the pasteur pipette (See Figure 3-3).




Figure 3-2





Figure 3-3

Cryotop® - Closed System

Thawing 1

Stand the CryotopSC on the Aluminum Block.




Thawing 2

Cut the marking point with Straw Cutter.

Put the cutting blade at the black marking point.

Turn the straw cap slowly to cut.




Thawing 3

Be ready to use pasteur pipette keeping the CryotopSC in liquid nitrogen. Set up the stop watch (with count up function). Check the time with the stop watch for the following steps.

Thawing 4

Insert the cut piece of the straw cap into the space between the CryotopSC and the Aluminum Block. This is to take out the CryotopSC easier.



Thawing 5 - For 1 minute

Quickly immerse the CryotopSC sheet into TS on the microscope stage by transferring it linearly. It should be within 1 second (See Figure 3-4). Find the Oocyte (Embryo) adjusting the focus on the black mark of the Cryotop sheet.

1 minute after immersing into TS, gently aspirate the Oocyte (Embryo) with the pasteur pipette after dispensing it from the sheet. Aspirate the Oocyte (Embryo) even if it does not dispense from the sheet. Also, aspirate TS until the Oocyte (Embryo) reaches 2mm from the tip of the pasteur pipette (See Figure 3-5).




Figure 3-4




Figure 3-5

THAWING



PART 4

DILUTION

Cryotop® and Cryotop®SC have different procedures.

Dilution - For 3 minutes

Blow out only **TS** in the pasteur pipette into the **BOTTOM** center of **DS** slowly (See Figure 4-1a), then gently place the Oocyte (Embryo) on the bottom of the **TS** layer (See Figure 4-1b). Leave it for 3 minutes. This is for mostly gradual displacement from **TS** to **DS**.




PART 5

PREPARATION FOR VITRIFICATION

Washing 1 - For 5 minutes

3 minutes later, after immersing into **DS**, gently aspirate the Oocyte (Embryo) in **DS** with the pasteur pipette. Also, aspirate **DS** until the Oocyte (Embryo) reaches 2mm from the tip of the pasteur pipette (See Figure 5-1).



Blow out only **DS** in the pasteur pipette into the **BOTTOM** center of **WS1** slowly (See Figure 5-2a), then gently place the Oocyte (Embryo) on the bottom there (See Figure 5-2b). Leave it for 5 minutes. This is also for mostly gradual displacement from **DS** to **WS1**.




Figure 5-2a




Figure 5-2b

Washing 2 - For 1 minute

5 minutes later, after immersing into **WS1**, aspirate the Oocyte (Embryo) with minimal volume of **WS1** with pasteur pipette (See Figure 5-3) and transfer it to the **TOP** center of **WS2**. After the Oocyte (Embryo) free-falls to the bottom of **WS2**, do the same work again in **WS2** (See Figure 5-4).




Figure 5-3





Figure 5-4


Washing 3

Transfer the Oocyte (Embryo) to a culture dish containing the appropriate culture medium. Incubate the Oocyte (Embryo) in a 37°C incubator to complete recovery. Completion of recovery: Oocyte (Embryo) for 2 hours for recommendation.



Vitrification Protocol

Embryo





KITAZATO

KITAZATO CORPORATION
100-10 Yanagishima, Fuji
Shizuoka 416-0932 Japan
Tel +81 545-65-7122
Fax +81 545-65-7128
contact@kitazato.co.jp

Dibimed

DIBIMED
C/ Jorge Comín, 3 Bajo 1-2
46015 Valencia, Spain
Tel +34 96 305 63 95
Fax +34 96 305 63 96
info@dibimed.com
Official Representative of Kitazato

