


Aspiration Pressure for different Needle Gauges for Kitazato OPU Needles



Needle Gauge Single Gauge	Aspiration Pressure	
	(mmHg)	(kPa)
16G	100 ± 25	13.3 ± 3.3
17G	125 ± 25	16.7 ± 3.3
18G	150 ± 25	20.0 ± 3.3
19G	175 ± 25	23.3 ± 3.3
20G	200 ± 50	26.7 ± 6.6
21G	225 ± 50	30.0 ± 6.6

Needle Gauge Double lumen	Aspiration Pressure	
	(mmHg)	(kPa)
16G	150 ± 25	13.3 ± 3.3
17G	175 ± 25	16.7 ± 3.3
18G	200 ± 25	20.0 ± 3.3