

**SUPPORTED AND
FULLY ECHOGENIC
ET CATHETER**

RYDEN

Experience the difference.
Clarity and control every time.



New

- **Fully echogenic** catheter with echo-tip on the guide
- **Outer metallic support** for more control
- **Bulb tip** for atraumatic transfer

Ultrasound image provided by Atsumi Yoshida, MD, PhD,
Kiba Park Clinic, JAPAN

KITAZATO®


RYDEN ET Catheter

Order Number	Description	Diameter (Fr)	Length (cm)
213312	Echogenic Supported Catheter with Guide	3	22
213355	Echogenic Supported Catheter with Guide	3	25
223312	Echogenic Supported Catheter and Guide with Obturator	3	22
223355	Echogenic Supported Catheter and Guide with Obturator	3	25

Order Number	Description
244017	Stylet for guide with 22cm catheter
244020	Stylet for guide with 25cm catheter



Catheter with metallic support that generates rigidity for a better control of the catheter



Echogenicity along the whole catheter, which will provide better visualization along the way to the uterine cavity



Bulb-tip with Echogenic ring that enhances a safe, smooth and atraumatic pass-way through the cervix



Pre-curved guide with memory to give an easy pass along



Large and user-friendly opening in the handle for insertion of the preloaded catheter into the guide