

Embryo Transfer Catheters

Designed following the advice of the most reputed gynaecologists



Straight ET Catheters

- Soft, flexible EC-PRO catheter with memory
- Echogenic EC-PRO Master model enhances visibility
- Supported catheter available
- Trial catheter for mock transfer procedure
- Stylet set with guide

Precurved ET Catheters

- Available in 3 and 4 French diameters
- Two tip options for gentle insertion: straight and bulb
- Sets include obturator
- Stylet available separately

Browse
our
catalogue



Straight ET


Type	Order Number	Description	Diameter (Fr)	Length (cm)
EC-PRO Routine	275018	Trial Catheter and Guide	4.7	18
	275023	Trial Catheter and Guide	4.7	23
	205118	Supported Catheter	4.7	18
	205123	Supported Catheter	4.7	23
	215018	Catheter and Guide	4.7	18
	215023	Catheter and Guide	4.7	23
	215118	Supported Catheter and Guide	4.7	18
	215123	Supported Catheter and Guide	4.7	23
	225523	Supported Catheter and Guide with Obturator	4.7	23



Type	Order Number	Description	Diameter (Fr)	Length (cm)
EC-PRO Master	205218	Echogenic Catheter	4.7	18
	205223	Echogenic Catheter	4.7	23
	215218	Echogenic Catheter and Guide	4.7	18
	215223	Echogenic Catheter and Guide	4.7	23
	205318	Echogenic Supported Catheter	4.7	18
	205323	Echogenic Supported Catheter	4.7	23
	215318	Echogenic Supported Catheter and Guide	4.7	18
	215323	Echogenic Supported Catheter and Guide	4.7	23



Type	Order Number	Description	Diameter (cm)	Quantity
Stylets	294015	Soft stylet with guide for 18cm catheter	15	10
	294020	Soft stylet with guide for 23cm catheter	20	10
	294115	Hard stylet with guide for 18cm catheter	15	10
	294120	Hard stylet with guide for 23cm catheter	20	10



Precurved ET



Type	Order Number	Description	Diameter (Fr)	Length (cm)
Bulb Tip Guide	213322	Catheter and Guide	3	22
	213325	Catheter and Guide	3	25
	214422	Catheter and Guide	4	22
	214425	Catheter and Guide	4	25
	223322	Catheter and Guide with Obturator	3	22
	223325	Catheter and Guide with Obturator	3	25
	224422	Catheter and Guide with Obturator	4	22
	224425	Catheter and Guide with Obturator	4	25
	223340	Catheter and Guide with Obturator	3	40



Type	Order Number	Description	Diameter (Fr)	Length (cm)
Straight Tip Guide	214322	Catheter and Guide	4	22
	214325	Catheter and Guide	4	25
	224322	Catheter and Guide with Obturator	4	22
	224325	Catheter and Guide with Obturator	4	25



Type	Order Number	Description	Diameter (Fr)	Length (cm)
Stylets	243017	Stylet for guide with 22cm catheter	3	17
	243020	Stylet for guide with 25cm catheter	3	20
	244017	Stylet for guide with 22cm catheter	4	17
	244020	Stylet for guide with 25cm catheter	4	20