

Oocyte Retrieval Needles



Main features

- Tip with smooth echogenic marking
- Thin wall design
- Sharp triple-cut blade
- Ergonomic handle with position marker
- Single and Double lumen alternatives

Browse
our
catalogue



LESS PAIN WITH ***THIN WALL DESIGN***

A Kitazato needle offers **the same inner diameter width for the aspiration** opposed to other brand thicker needles.

This allows oocyte retrieval to be carried out with a smaller gauge needle and therefore reduces the possibility of bleeding and the risk of injuries to the uterine lining, whilst minimizing stress in the oocyte during aspiration.



Kitazato Needle
Less damage to the ovary



Competitor Needle

PRODUCT REFERENCES

Type	Order Number	Gauge	Length (cm)	Aspiration Line (cm)*
Single Lumen Pump Aspiration	326350	16	35	80
	327350	17	35	80
	328350	18	35	80
	329350	19	35	80
	320350	20	35	80
	320300	20	30	80
	321300	21	30	80

Type	Order Number	Gauge	Length (cm)	Aspiration Line (cm)*
Single Lumen Pump Aspiration with Vacuum Line	336350	16	35	80 + 60
	337350	17	35	80 + 60
	338350	18	35	80 + 60
	339350	19	35	80 + 60
	330300	20	30	80 + 60

*Other aspiration line length available

Type	Order Number	Gauge	Length (cm)	Aspiration Line (cm)
Single Lumen with Luer Lock	316001	16	32.5	-
	317001	17	32.5	-
	318001	18	32.5	-
	319001	19	32.5	-
	310001	20	32.5	-

Type	Order Number	Gauge	Length (cm)	Aspiration Line (cm)	Flushing Line (cm)
Double Lumen	376350	16	35	80	100
	377350	17	35	80	100
	378350	18	35	80	100

Type	Order Number	Product Name	Gauge	Length (cm)	Aspiration Line (cm)
Special Procedure	268003	Towako Method	18	32.5	-
	367350	Manual Aspiration Set for Syringe	17	35	60
	357353	Aspiration Set with Stopcock	17	35	80 + 40 with stopcock
	347275	IVM Needle	17 & 20	27.5 & 35	80