

# NEW IVF MEDIA



## Gamete Buffer

- For flushing and retrieval of oocytes
- To be used in ambient atmosphere for gamete handling and manipulation
- Ready to use
- It can be used up to 7 days after opening

Type	Order Number	Code	Contents
Gamete Buffer	95720	GBSD-20x5	5x20mL
	95721	GBSD-50	50mL
	95722	GBSD-100	100mL

### COMPONENTS

- Calcium Chloride
- Gentamicin
- Glucose
- HEPES
- Human Serum Albumin (4.0g/l)
- Magnesium Sulfate
- Potassium Chloride
- Sodium Phosphate
- Sodium Bicarbonate
- Sodium Chloride
- Lactate
- Pyruvate

### QUALITY CONTROL

pH 7.3-7.6 / Osmolality 270-290mOsm/kg  
 Endotoxin <0.25EU/mL / Sterility Test  
 Mouse Embryo Assay ≥ 80%  
 Storage: 2-8 °C  
 Shelf Life: 18 months

