

Intrauterine Insemination Catheters



Main features

- Double lateral opening tip
- Smooth and rounded tip
- Four types of rigidity
- Integrated inner stylet

Browse
our
catalogue



Type	Order Number	Description	Diameter (Fr)	Length (cm)
IUI	131194	Soft	5.2	18
	131192	Intermediate	6	18
	131196	Ultimate (with built-in stylet)	6	18
	131191	Hard (with built-in stylet)	6	18

SOFT IUI CATHETER

This soft catheter minimizes trauma to the patient.



INTERMEDIATE IUI CATHETER

Allows for some flexibility and malleability.



ULTIMATE IUI CATHETER

Malleable catheter that adopts the desired shape to facilitate its insertion.

HARD IUI CATHETER

Rigid catheter with built-in stylet. Adopts and maintains the desired shape.

