

SUPPORTED CATHETER WITH ECHO-TIP

RYDEN™

Experience the difference.
Clarity and control every time.



New

- Catheter with **echo-tip on the guide**
- **Outer metallic support** for more control
- **Bulb tip** for atraumatic transfer

KITAZATO®

PRODUCT REFERENCES

RYDEN™ ET Catheter

Order Number	Description	Diameter (Fr)	Length (cm)
213342	Supported Catheter and Guide	3	22
213365	Supported Catheter and Guide	3	25
223342	Supported Catheter and Guide with Obturator	3	22
223365	Supported Catheter and Guide with Obturator	3	25

Order Number	Description
244027	Stylet for guide with 22cm catheter
244030	Stylet for guide with 25cm catheter



Catheter with metallic support that generates rigidity for a better control of the catheter



Bulb-tip with Echogenic ring which will provide better visualization while enhancing comfort and safety for the patient.



Pre-curved guide with memory to give an easy pass along



Large and user-friendly opening in the handle for insertion of the preloaded catheter into the guide