

NEW IVF MEDIA



Scan for
more info

Sperm Freeze

- **Ready-to-use media intended** for cryopreservation of human sperm.
- For samples with **good initial quality**.
- **Extended Stability:** Maintains efficacy for 7 days after opening.

Order Number	Code	Content
95951	SFD2-10	10ml
95962	SFD2-20x5	5x20ml

COMPONENTS

- HEPES buffer
- Physiological salts
- Glycine
- Glucose
- Lactate
- Glycerol (15%)
- Sucrose
- Human Serum Albumin (4.0 g/l)

QUALITY CONTROL

pH 7.20-7.60
Sterility Test
Endotoxins: <0.25 EU/ ml
Human Sperm Survival Assay ≥80%
Storage: 2-8°C
Shelf Life: 18 months

KITAZATO®



Scan for
more info

Sperm Freeze G-Plus

- **Ready-to-use media used** for cryopreservation of human sperm.
- **Enhanced Glycerol Content (27%)**: for compromised samples.
- **Ideal for long-term** storage and sperm banks.

Order Number	Code	Content
95971	SFD4-10	10ml
95980	SFD4-1x5	5x1ml

COMPONENTS

- HEPES buffer
- Physiological salts
- Glycine
- Glucose
- Lactate
- Glycerol (27%)
- Sucrose
- Human Serum Albumin (4.0 g/l)

QUALITY CONTROL

pH 7.20-7.60
Sterility Test
Endotoxins: <0.25 EU/ ml
Human Sperm Survival Assay ≥80%
Storage: 2-8°C
Shelf Life: 18 months

KITAZATO CORPORATION
100-10 Yanagishima, Fuji
Shizuoka 416-0932 Japan
Tel +81 545-65-7122
Fax +81 545-65-7128
contact@kitazato.co.jp

KITAZATO | DIBIMED
C/ Jorge Comin, 3 Bajo 1-2
46015 Valencia, Spain
Tel +34 96 305 63 95
Fax +34 96 305 63 96
info@dibimed.com
Official Representative of Kitazato

KITAZATO®